

## Introduction

Excellent optical system, high level mechanical system, advanced circuit control system, rigorous production process, friendly and intuitive software interface, good technical specifications, stable and reliable performance can meet the analysis requirements from high level and professional customers.



## Main Features

### ***Appearance and internal structure***

Modern and elegant appearance, extendable design, separate structure design for optical and circuit system can efficiently avoid the loss of photometric energy.

### ***Convenient and intuitive operation interface***

This series has 7-inch high resolution colored capacitive touch screen and newly developed UV-SUPER2.0 software with strong functions, which make the operation simple and easy.

### ***Excellent performance and stability***

Totally enclosed monochromator and optical mirror coated with SiO<sub>2</sub> guarantee the optical components are not influenced by gas and environment.

- ① Philips and Hamamatsu lamps.
- ② Newly improved screw pole drive structure makes good wavelength repeatability and high wavelength accuracy.
- ③ Totally new design, superior materials and rigorous production process.

### ***Advanced photoelectric test system***

- ① 32 bit ARM11 microcontroller with clock speed up to 533MHz.
- ② 20 bit analog digital device specialized for photoelectric data collection and processing from BB company.
- ③ Support internal huge data storage, mouse operation and big SD card memory.

### ***Simple and convenient maintenance***

- ① Socket type lamps make the optical adjustment not necessary and maintenance much easier.
- ② Separated Optical and circuit system has no cross influence and make the instrument more reliable.

# T-9100/9200

## Specifications

| MODEL                          | T-9100  | T-9200                 | T-9200S                | T-9200A                |
|--------------------------------|---|------------------------|------------------------|------------------------|
| Display                        | 7 inch TFT colored capacitive touch screen        |                        |                        |                        |
| Wavelength Range               | 190 - 1100nm                                      |                        |                        |                        |
| Optical System                 | Single Beam                                       | Double Beam            |                        |                        |
| Spectral Bandwidth             | 2nm   | 2nm                    | 1nm                    | 0.5,1,2,4nm            |
| Wavelength Accuracy            | ±0.5nm  | ±0.5nm                 | ±0.3nm                 | ±0.3nm                 |
| Wavelength Repeatability       | ≤0.2nm  | ≤0.2nm                 | ≤0.1nm                 | ≤0.1nm                 |
| Photometric Range              | 0-200%T, -0.3-3.0A, 0-9999C (0-9999F)             |                        |                        |                        |
| Photometric Accuracy           | 0.2%T (0-100%T), ±0.002A(0-0.5A), ±0.004A(0.5-1A) |                        |                        |                        |
| Photometric Repeatability      | ≤0.15%T (0-100%T), 0.001A(0-0.5A), 0.002A(0.5-1A) |                        |                        |                        |
| Scanning Speed                 | Low, Medium, High (up to 3000nm/min)              |                        |                        |                        |
| Stray Light                    | ≤0.05%T@220nm,360nm                               |                        |                        |                        |
| Baseline Flatness              | ±0.003A   | ±0.002A                | ±0.001A                | ±0.001A                |
| Drift                          | 0.003A/30min<br>@500nm                            | 0.002A/30min<br>@500nm | 0.001A/30min<br>@500nm | 0.001A/30min<br>@500nm |
| Noise                          | 0.0005A@500nm                                     |                        |                        |                        |
| Working Mode                   | T,A,C,E   |                        |                        |                        |
| Wavelength Setting             | Automatic   |                        |                        |                        |
| Detector                       | Solid Silicon Photodiode                          |                        |                        |                        |
| Light Source                   | Tungsten Halogen/Deuterium Lamp                   |                        |                        |                        |
| Output Port                    | USB HOST, USB DRIVE, RS232                        |                        |                        |                        |
| Power Requirements             | AC 110-220V 50-60Hz                               |                        |                        |                        |
| Humidity Range                 | Less Than 85%                                     |                        |                        |                        |
| Shipping Dimensions and Weight | 770*630*340mm, 27kg                               |                        |                        | 880*690*530mm,<br>45kg |



Absorbency and transmittance test

## Photometry



To test sample solution concentration, you can choose different methods like coefficient, standard curve, linearity, linear zero crossing and quadratic. Operators can choose single, double and tri-wavelength and change the coefficients of double and tri-wavelength. Advanced arithmetic makes curvilinear regression more precise and test data more accurate.

## Quantitative Measurement



To test the sample chemical reaction process by fixed time scanning the sample solution with fixed wavelength. The equipment can calculate its changing rate after entering the corresponding parameters.

## Kinetics Measurement(Time Scanning)



To test sample solution absorbency peak, can scan the sample characteristic curve of any wavelength range between 190 and 1100nm. And do the chart overlay and arithmetic.

## Wavelength Scanning(Qualitative Test)

| Options | 190   | 210   | 300   | 500   | 700   | 800   | 900   | 1000  | 1100 (nm) |
|---------|-------|-------|-------|-------|-------|-------|-------|-------|-----------|
| 01      | 0.210 | 0.410 | 0.210 | 0.210 | 0.280 | 0.910 | 0.210 | 0.210 | 0.210     |
| 02      | 0.210 | 0.410 | 0.210 | 0.210 | 0.280 | 0.910 | 0.210 | 0.210 | 0.210     |
| 03      | 0.210 | 0.410 | 0.210 | 0.210 | 0.280 | 0.910 | 0.210 | 0.210 | 0.210     |

It is much more convenient for users to test the absorbency or do the arithmetic in case of several wavelengths for the same sample solution, which is much simpler than single wavelength testing.

## Multi Wavelength Measurement

| Options | 1100 Abs | 600 Abs | 400 Abs | 比率    | 蛋白质(浓度) | 核酸(浓度) |
|---------|----------|---------|---------|-------|---------|--------|
| 01      | 0.210    | 0.410   | 0.210   | 0.910 | 0.210   | 0.210  |
| 02      | 0.210    | 0.410   | 0.210   | 0.910 | 0.210   | 0.210  |
| 03      | 0.210    | 0.410   | 0.210   | 0.910 | 0.210   | 0.210  |

It is a special function for specific users and make the operation easier.

## DNA/Protein Measurement

# Accessories



Manual 4-position cell holder (standard for single beam)



Single-hole cuvette holder (Standard for Double Beam)



Automatic 8-position round cell holder



Manual 4-position 10cm cell holder



Manual 4-position film holder



Single hole film holder



Single hole long optical path holder



Tube rack



Adjustable XY micro cell holder

# Comparison Table

|          | UV/Vis.      | Optical System | Display      | Slit Width   | Wavelength Accuracy | Wavelength Repeatability | Stray Light          | Light Source                     | Page  |         |             |      |      |                                  |             |      |      |              |      |              |        |
|----------|--------------|----------------|--------------|--------------|---------------------|--------------------------|----------------------|----------------------------------|-------|---------|-------------|------|------|----------------------------------|-------------|------|------|--------------|------|--------------|--------|
| E-1000V  | Vis.         | Single         | 70*40mm LCD  | 4 nm         | ±2 nm               | ≤1 nm                    | ≤0.15%T@360nm        | Tungsten Halogen Lamp            | 1/2   |         |             |      |      |                                  |             |      |      |              |      |              |        |
| E-1000UV | UV           |                |              |              |                     |                          |                      | Tungsten Halogen /Deuterium Lamp |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| C-7000V  | Vis.         | Single         | 7 inch FTF   | 2 nm         | ±0.3 nm             | ≤0.1 nm                  | ≤0.05%T @220nm,360nm | Tungsten Halogen Lamp            | 3/4   |         |             |      |      |                                  |             |      |      |              |      |              |        |
| C-7000UV | UV           |                |              |              |                     |                          |                      | Double                           |       | 1 nm    | 0.5,1,2,4nm | 2 nm | 1 nm | Tungsten Halogen /Deuterium Lamp |             |      |      |              |      |              |        |
| C-7100   |              |                |              |              |                     |                          |                      |                                  | 1 nm  |         |             |      |      |                                  | 0.5,1,2,4nm | 2 nm | 1 nm | 0.5,1,2,4 nm |      |              |        |
| C-7100S  |              |                |              |              |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              | 1 nm | 0.5,1,2,4 nm | 1, 2nm |
| C-7100A  |              |                |              |              |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| C-7200   | 1 nm         | 0.5,1,2,4 nm   | 1, 2nm       |              |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| C-7200S  | 0.5,1,2,4 nm | 1, 2nm         |              |              |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| C-7200A  | 1, 2nm       |                |              |              |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| C-7200PC |              |                |              |              |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| X-8200   | UV           | Double         | 7 inch FTF   | 2 nm         | ±0.6 nm             | ≤0.2 nm                  | ≤0.1%T @220nm,360nm  | Xenon Lamp                       | 9/10  |         |             |      |      |                                  |             |      |      |              |      |              |        |
| X-8200S  |              |                |              | 1 nm         |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| X-8200A  |              |                |              | 0.5,1,2,4 nm |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| X-8200T  |              |                | Touch Screen | 2 nm         |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| X-8200TS |              |                |              | 1 nm         |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| X-8200TA |              |                |              | 0.5,1,2,4nm  |                     |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| T-9100   | UV           | Single         | Touch Screen | 2 nm         | ±0.5 nm             | ≤0.2 nm                  | ≤0.05%T @220nm,360nm | Tungsten Halogen /Deuterium Lamp | 11/14 |         |             |      |      |                                  |             |      |      |              |      |              |        |
| T-9200   |              | Double         |              | 1 nm         |                     |                          |                      |                                  |       | ±0.3 nm | ≤0.1 nm     |      |      |                                  |             |      |      |              |      |              |        |
| T-9200S  |              |                |              |              | 0.5,1,2,4,nm        |                          |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |
| T-9200A  |              |                |              |              |                     | 0.5,1,2,4,nm             |                      |                                  |       |         |             |      |      |                                  |             |      |      |              |      |              |        |